

pipeline

December 24, 2020

1 BCI-4-ALS

1.1 Mid-Semester Project

1.1.1 Project Overview

Over the mid-semester project we have been integrated the different parts of BCI which been discussed in the course. The project focused on Motor Imagery interface which aims to classify between imagery of moving right hand, left hand and idle.

Over the first part in the report we present the different codes we made for the MI classification. This part can be used as a README file for next years.

```
[1]: import yaml
import pyxdf
import os
import numpy as np
import mne
import pickle
```

Configurations All the scripts using the `config.yaml` configuration file. The first part of each script include loading the config file and extract the parameters which relevant for the current part. The configuration file arranged to be easy to read and change. For example, we can get and set LabStreamingLayer parameters by:

```
[2]: config = yaml.full_load(open('../MI/config.yaml', 'r'))
lsl_params = config['lsl']

print('\n'.join('{}: {}'.format(k, v) for k, v in lsl_params.items()))
```

```
start_experiment: 1111
start_trial: 111
end_trial: 999
end_experiment: 9999
```

The configuration file centralize all the interface parameters into one file.

MI1 - Record Experiment The `MI1_record_experiment.py` script responsible for the visualization of the experiment and recording it. The main function `record_experiment()` get paradigm

(for now MI or SSVEP) and running experiment according given paradigm. For now we only partially implemented the SSVEP paradigm and focused on MI. We also implemented an interactive GUI which help the user for running an experiment.

The `MI_record()` function is the main function for recording a MI experiment. The function create folder for the subject, initialize LSL stream, creating random stimulus vector and start presenting the subject the images. The function push before and after each image a marker to indicate the trial start and end. Parameters for this part can be change in the config file, such as the amount of trials, length of each trial, etc.

The output of MI1 is xdf file which contains the EEG and markers record. The second file is the stimulus vector which contains the order of the images the subject saw in the experiment (i.e., the y vector). The files' names are EEG.xdf and stimulus_vector.csv and both of the save in the subject directory.

MI2 - Data Pre-Processing The `MI2_preprocess.py` script responsible for getting the subject's folder directory from the config file and for each day of recording cleans the data. Currently, the pre-processing script using low-pass, high-pass and notch filters to clean the data. The parameters for the filters are also part of the config file under `config['preprocess']['filter']`.

The output of MI2 script is `EEG_clean.csv` file, which contains the EEG data after the pre-processed. Additionally, the script export a `.info` json file. The `.info` file contains all the info about the EEG stream. MI4 uses the info of the streaming. The output is automatically saves in the data directory.

A plot of the clean EEG data seem as follows:

```
[8]: #----- todo: find out if this is a normal psd plot -----

%matplotlib inline
# Params
eeg_clean_path = '../data/noam/2/EEG_clean.csv'
ch_names = ['C03', 'C04', 'P07', 'P089', 'O01', 'O02', 'F07', 'F08', 'F03', 'P04', 'T07', 'T08', 'P03']
s_rate = 125

# Create mne info & raw
eeg_clean = np.genfromtxt(eeg_clean_path, delimiter=',', skip_header=1)[: , 1:]
info = mne.create_info(ch_names, s_rate)
raw_eeg_clean = mne.io.RawArray(eeg_clean.T, info)

# Plot
fig = raw_eeg_clean.plot_psd(picks=ch_names, show=False)
```

Creating RawArray with float64 data, n_channels=13, n_times=95705

Range : 0 ... 95704 = 0.000 ... 765.632 secs


Ready.

Effective window size : 16.384 (s)

<ipython-input-8-f08db6f710a5>:15: RuntimeWarning: Channel locations not

available. Disabling spatial colors.

```
fig = raw_eeg_clean.plot_psd(picks=ch_names, show=False)
```



MI3 - Data Segmentation In the `MI3_segment_data.py` script we split the data for each trial. The start and end of each trial was according to the markers streaming we created while recording the EEG. The script get the subject folder from the config file (under `config['data']['subject_folder']`) and split for trial each EEG record in each day.

The output of MI3 script is a pickle file named `EEG_trials.pickle` which located in the corresponded day directory. The file is a list with ndarray where each ndarray is the corresponded EEG data of the trial.

```
[5]: eeg_trials_path = '../data/noam/3/EEG_trials.pickle'
eeg_trials = pickle.load(open(eeg_trials_path, 'rb'))

print('Number of trials: {} \n Trials dimensions: {}'.format(len(eeg_trials),
    ↳ eeg_trials[0].shape))
```

Number of trials: 120

Trials dimensions: (616, 13)

And a PSD plot of various trial look like:

```
[16]: eeg_trial = mne.io.RawArray(eeg_trials[103].T, info)
fig = eeg_trial.plot_psd(picks=ch_names)
```

Creating RawArray with float64 data, n_channels=13, n_times=614


Range : 0 ... 613 = 0.000 ... 4.904 secs

Ready.

Effective window size : 4.912 (s)

<ipython-input-16-6dbe10178dc0>:2: RuntimeWarning: Channel locations not available. Disabling spatial colors.

```
fig = eeg_trial.plot_psd(picks=ch_names)
```



MI4 - Feature Extraction